

All
Clinicians

Epic

Image Viewer




Overview of Feature/Changes

Effective Wednesday March 13, 2019 at 7:00AM, clinicians will have the ability to view images without having to leave hyperspace. This can be done using a new viewer. You can choose to view one image or all images for a specific patient. The viewer is available in the Ambulatory, Inpatient and Emergency Department context.

Accessing the Viewer

Ambulatory, Inpatient, Emergency Department: Two options

1. From Chart Review, Click on the Imaging tab.



- Choose the Imaging study you would like to view.
- Click on the blue hyperlink in the report pane under PACS Images.

02/19/2019		ELECTRORETINOGRAPHY (ERG) -...	
12/12/2018	ULT379956632	US LIMITED (IN OFFICE)	
12/07/2018	1032139	INTERVENTIONAL RA... IR IMAGE GUIDED NEEDLE DRAIN...	
11/20/2018		DEXA - VERTEBRAL FRACTURE AS...	
11/20/2018		DEXA HIP AND OR SPINE 77080	
11/20/2018		DEXA TOTAL BODY 77080	
11/20/2018		DEXA WB/ST 77081	

PACS Images
[Show images for US Limited \(in office\)](#)


- This will open the viewer in hyperspace for that specific study.



2. In Chart Review when in the Imaging tab, **Click** the Imaging button.




- This will open the viewer in hyperspace and allow the clinician to view all imaging studies associated for the patient.



Inpatient/Emergency Department: More Activities Option

In the Inpatient context, clinicians can also access SynMobility through the Imaging activity. To do this, **Click** on the More button and **Click** on Imaging. You can make this a favorite by **Clicking** the Star. You can adjust where to place it in the activity tabs using the customize feature.



Emergency Department: Provider Trackboard Option

In the Emergency department context, clinicians can also access the Imaging button from the Provider Trackboard.

- **Click** on the Workup tab.
- In the Imaging section, **Click** View Image.

My Patients | Waiting for Provider | All Patients (11) | Waiting Room | Consults

My + Unassigned | Rad Tracking | To Be Admitted | Alerts | Care in Triage

Results Pending | A | B | Finish Up Fast

Bed	Patient	Age	Comp	ESI	TT	RN	MD	RES	New I	Img C	Img C	Requ	Bed S	Comr
B18	Expect, Uh (27 y.o.) (F)	27...			...	S...	A...							
H0	Francisalbert, Sinatra (89 y...	89...	sick		...	T...								
B20	Lab, Validation (30 y.o.) (M)	30...	te...		...	T...	C...							
B19	Legal, Status (30 y.o.) (M)	30...	va...		...	D...								
A07	Synapse, Mobility (24 y.o.)...	24...	A...		...	T...								
A33	Test, Admit (16 y.o.) (M)	16...			...	T...	S...						C...	R...
2R	Test, Cimuhed (15 y.o.) (F)	15...		2	...	T...	Jo...							
A11	Test, Peduhsep (8 y.o.) (F)	8...			...	T...								
3R	Test, Trauma (61 y.o.) (F)	61...			...	T...								
H37	Testgc, Rjdnew (30 y.o.) (M)	30...			...	T...								
A12	Today, Legal (35 y.o.) (M)	35...			...	T...								

BestPractice Advisories

First Provider Eval

01/30 0806

ED First Provider Evaluation

First Provider Eval FILE

Labs

ECG

Imaging **NEW**

XR Humerus AP Lateral Left

Image Available 01/30/19 0925

XR Shoulder Complete Left

Final Result 01/30/19 0934

View Image Interpret

View Image

Accessing images from the Epic toolbar

You can also access images in Synapse by using the Synapse Imaging button on the Epic toolbar.

1. In the Epic toolbar, Click Synapse Imaging.

Hyperspace - MEDICINE - POC - ATTENDING I.

60 - Hospital Chart Completion

Epic

Patient Lists | In Basket | Patient Lookup | Schedule | Prep for Surg/Procedure | Orders for Admission | UpToDate | SmartList | Personalize | Carefusion | Remind Me | **Synapse Imaging**

Patient Lists

2. This will open the Synapse Mobility application window. Enter your Username and Password.
 - This avenue will allow you to search Syntemp images or search additional criteria for non-linked patients using MRN numbers.

5/7/2019 Today

https://synmobility.upstate.edu:8443/ Synapse Mobility Login

SYNAPSE

SYNAPSE mobility

Enter Name:

Enter Password: