

# High Sensitivity Cardiac Troponin T (hs-cTnT)



How to use and comparison to existing troponin  
testing


# Why do we need High sensitivity cardiac Troponin (hs-cTnT)?

- Current troponin assays cannot accurately measure low levels of troponin – we cannot detect trends at lower levels or increases warning of cardiac ischemia
- High sensitivity Troponin can accurately measure low levels of troponin

## Comparison of measuring range

Current cTnT (or 4<sup>th</sup> Gen)

hs-cTnT (or 5<sup>th</sup> Gen)



Clinical Biochemistry 2015;48:201-3

<https://www.aacc.org/publications/clin/clin-stat/2015/april/16/New-Guidance-on-High-sensitivity-Cardiac-Troponin-Testing>

## This increased sensitivity allows us to:

- Diagnose early MI
- Reduce the time to rule-in and rule-out acute MI
- Improve clinical accuracy of testing
- Improve MI patient outcomes
- Reduce turnaround time in the ED for MI rule-out patients

Troponin T High Sensitive

Process Instructions: Tier 1 (all credentialed providers)

Frequency: ONCE

Starting: 6/1/2021 Today Tomorrow At: 1525

First Occurrence: Today 1525

Show Scheduled Times

06/01/21 1525

Reference Links: Tiered Categorization

Comments: Add Comments (F6)

Next Required Link Order

# hs-cTnT provides greater sensitivity for MI in women

- Current Troponin assay has a single cut-off for both sexes (10ng/L) for suspicion of MI
- The increased sensitivity of hs-cTnT allows use of a lower cutoff for female patients
- Using hs-cTnT we can use sex-specific thresholds for risk of cardiac damage
  - Normal female patients have less than 14 ng/L of Troponin as measured by hs-cTnT
  - Normal male patients have less than 22 ng/L of Troponin as measured by hs-cTnT
- To ensure that no MIs are missed, we use thresholds of 10 ng/L (females) and 15 ng/L (males) to rule out MI on the initial testing
- As with current Troponin assay, patients are best evaluated with serial hs-cTnT testing

# How to use hs-cTnT testing – Acute MI

For Rule Out/Rule In Patients Suspicious for Acute MI

Time 0 hs-cTnT Epic Code: LAB9290

< 10 ng/L Female  
< 15 ng/L Male  
AMI ruled out

> 10 ng/L Female  
> 15 ng/L Male  
Draw 2<sup>nd</sup> hs-cTnT at 2 hours from time 0

Increase in hs-cTnT < 3 ng/L  
AMI ruled out

Increase in hs-cTnT > 3 ng/L and < 10 ng/L  
Draw 3<sup>rd</sup> hs-cTnT at 6 hours from time 0

Increase in hs-cTnT > 10 ng/L  
AMI ruled in

Increase in hs-cTnT < 12 ng/L from time 0  
AMI ruled out

Increase in hs-cTnT > 12 ng/L from time 0  
AMI ruled in


○ All rule-ins and rule-outs should be further stratified with history, EKG, age, and other risk factors

○ **Critical value** (called to clinical team by lab) = **100 ng/L**

# How to use hs-cTnT values – Chronic MI


- Distinguish **acute** from **chronic** MI → look for significant change (rise or fall) of hs-cTnT results
  - ❖ Troponin level changes from the time of ischemia and type of injury (see graph)\*
  - ❖ “Significant change” (rise and/or fall) definition depends if initial hs-cTnT is >10 ng/L (females) or >15 ng/L (males)
    - If the initial result is **below** threshold: a change of at least **50%** is significant
    - If the initial result is **above** threshold: a change of at least **20%** is significant

- Collect 2<sup>nd</sup> sample 2 hrs after 1<sup>st</sup> sample to determine the pattern of hs-cTnT change (and additional samples if indicated)
- Always interpret hs-cTnT results along with clinical picture (History, ECG, HEART Score, etc.)



# What are clinical conditions that can cause chronic hs-cTnT elevations ?

- Ischemic or non-ischemic heart failure patients with different forms of cardiomyopathy
- Myocarditis,
- Heart contusion,
- Pulmonary embolism
- Stroke
- Subarachnoid hemorrhage
- Hypertensive crisis
- Renal failure →
- Sepsis
- Diabetes
- Drug-induced cardiotoxicity
- Critical illness



## Differential Dx between **acute** and **chronic hs-cTnT** elevations,

- A serial sampling to observe a rise/fall of hs-cTnT above threshold (14 ng/L for females, 22 ng/L for males) **AND**
- Consistent with the clinical assessment, including ischemic symptoms and electrocardiographic changes

# hs-cTnT vs cTnI on iSTAT

- The point-of-care iSTAT device cardiac troponin test is a **conventional** troponin test.
- The testing results between the iSTAT and the laboratory hs-cTnT are **not** equivalent and should **not** be used interchangeably or as part of the trending to determine acute vs chronic cardiac damage.
- iSTAT POC troponin testing should **only** be used when there is not time to obtain the hs-cTnT testing from the core laboratory.
- If a patient result(s) of iSTAT disagree with result of hs-cTnT tested at core-laboratory, the hs-cTnT result should be used for patient care interpreted with other signs/symptoms.

Apple et al. 2021; Christenson et al. 2017



**If testing results do not match patient presentation or for any other questions,  
contact the Chemistry Laboratory directors for assistance at: 315-464-9175**

Email: [caoz@upstate.edu](mailto:caoz@upstate.edu) or  
[elkinsm@upstate.edu](mailto:elkinsm@upstate.edu)

# Select References

- Apple FS, Fantz CR, Collinson PO 2021 Implementation of High-Sensitivity and Point-of Care Cardiac Troponin Assays into Practice: Some Different Thoughts. *Clin Chem* 67(1): 70-8
- Christenson RH, Jacobs E, Uettwiller-Geiger D, Estey MP, Lewandrowski K, Koshy TI, Kupfer K, Li Y, Wesenberg JC 2017 Comparison of 13 Commercially Available Cardiac Troponin Assays in a Multicenter North American Study. *JALM* 1(5): 544-61
- Boeddinghaus J, Nestelberger T, Twerenbold R, Koechlin L, Meier M, Troester V, Wussler D, Badertscher P, Wildi K, Puelacher C, du Fay de Lavallaz J, Gimenez MR, Zimmermann T, Hafner B, Potlukova E, Miro O, Martin-Sanzchez FJ, Keller DI, Reichlin T, Mueller C, APACE Investigators 2019 High-sensitivity Cardiac Troponin I Assay for Early Diagnosis of Acute Myocardial Infarction. *Clin Chem* 65(7); 893-904
- Giannitsis E, Kurz K, Hallermayer K, Jarausch J, Jaffe AS, Katus HA 2010 Analytical validation of a high-sensitivity cardiac troponin T assay. *Clin Chem* 56(2); 254-61
- Kraus D, von Jeinsen B, Tzikas S, Palapies L, Zeller T, Bickel C, Fette G, Lackner KJ, Drechsler C, Neumann JT, Baldus S, Blankenberg S, Munzel T, Wanner C, Zeiher AM, Keller T. Cardiac Troponins for the Diagnosis of Acute Myocardial Infarction in Chronic Kidney Disease. *J Am Heart Assoc* 7(19): Epub
- Neumann JT, Twerenbold R, Ojeda F, Sorensen NA, Chapman AR, Shah ASV, Anand A, Boeddinghaus J, Nestelberger T, Badertscher P, Mokhtari A, Pickering JW, Troughton RW, Greenslade J, Parsonage W, Mueller-Hennesen M, Gori T, Jernberg T, Morris N, Liebetrau C, Hamm C, Katus HA, Munzel T, Landmesser U, Salomaa V, Iacoviello L, Ferrario MM, Giampaoli S, Kee F, Thorand B, Peters A, Borchini R, Jorgensen T, Soderberg S, Sans S, Tunstall-Pedoe H, Kuulasmaa K, Renne T, Lackner KJ, Worster A, Body R, Ekelund U, Kavsak PA, Keller T, Lindahl B, Wild P, Giannitsis E, Than M, Cullen LA, Mills NL, Mueller C, Zeller T, Westermann D, Blakenberg S. Application of High-Sensitivity Troponin in Suspected Myocardial Infarction. *N Engl J Med* 380(26): 2529-40
- Reichlin T, Cullen L, Parsonage WA, Greenslade J, Twerenbold R, Moehring B, Wildi K, Mueller S, Zellweger C, Mosimann T, Gimenez MR, Rentsch K, Osswald S, Muller C 2015 Two-hour algorithm for triage toward rule-out and rule-in of acute myocardial infarction using high-sensitivity cardiac troponin T. *Am J Med* 128(4); 369-79
- Reidlinger D, Mockel M, Muller C, Holert F, Searle J, von Recum J, Slagman A 2018 High-Sensitivity cardiac troponin T for diagnosis of NSTEMI in the elderly emergency department patient: a clinical cohort study. *Biomarkers* 23(^); 551-7
- Thygesen K, et al. writing for ESC/ACF/AHA/WHF Task Force 2012 Third Universal Definition of Myocardial Infarction. *JACC* 60(16): 1581-98
- Von Jeinsen B, Brandebussemeyer S, Grued D, Bormann J, Psyraakis D, Wolter JS, Kriechbaum S, Doerr O, Troidl C, Nef H, Hamm CW, Liebetrau C, Keller T 2020 Evolution of the effect of older age on the diagnostic specificity of high-sensitivity troponin I and high-sensitivity troponin T in patients with suspected acute myocardial infarction. *ESC Congress 2020 Abstract*