

Aldosterone, Renin, CMV IgM testing change

- Starting **September 10th 2026**, testing for Aldosterone (**LAB123152**), Renin (**LAB123153**), and Renin:Aldosterone ratio (**LAB123154**), and CMV IgM (**LAB123155**) will transition from send-out testing to in-house testing.
- Performing this testing in-house will improve turn-around time and improve consistency
- The normal ranges for the in-house assays differ from the send-out assay

Reference intervals change:

	Old	New
Aldosterone	0.0-30.0 ng/dL	Upright <39.3 ng/dL; Supine <23.6 ng/dL
Renin activity	0.167-5.38 ng/mL/hr	Upright/supine and >/≤40 yrs normal ranges, pg/mL
Aldosterone/Renin ratio	0.0-20.0	< 3.7
CMV IgM antibodies	POS/NEG	<30.0 AU/mL negative; ≥30 AU/mL and <35 AU/mL equivocal; ≥ 35 AU/mL positive

- For further guidance on interpretation and usage of this testing, please see:
 SUNY Upstate Laboratory Test Catalog: <https://www.upstate.edu/pathology/healthcare/orders.php?search=>
 Primary Aldosteronism: An Endocrine Society Clinical Practice Guideline 2025, <https://doi.org/10.1210/clinem/dgaf284>

For any questions, please contact the Core Laboratory (315-464-4459)